

For Research Use Only

POLI Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60889-1-PBS



Basic Information

Catalog Number: 60889-1-PBS	GenBank Accession Number: BC032662	Purification Method: Protein G purification
Concentration: 1 mg/ml	GeneID (NCBI): 11201	CloneNo.: 2B2G1
Source: Mouse	UNIPROT ID: Q9UNA4	
Isotype: IgG1	Full Name: polymerase (DNA directed) iota	
Immunogen Catalog Number: AG4563	Calculated MW: 715 aa, 80 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51303-1, POLI Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60889-1-PBS. Detection antibody: 60889-2-PBS. Standard:Ag4563. Range: 0.781-200 ng/mL.