

For Research Use Only

# MASTL Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60891-1-PBS



## Basic Information

Catalog Number:

60891-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG9662

GenBank Accession Number:

BC009107

GeneID (NCBI):

84930

UNIPROT ID:

Q96GX5

Full Name:

microtubule associated  
serine/threonine kinase-like

Calculated MW:

879 aa, 97 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1D3H3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP51306-1, MASTL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60891-1-PBS. Detection antibody: 60891-2-PBS. Standard:Ag9662. Range: 0.781-100 ng/mL.