

For Research Use Only

Nectin 3 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60892-2-PBS



Basic Information

Catalog Number:

60892-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28680

GenBank Accession Number:

BC017572

GeneID (NCBI):

25945

UNIPROT ID:

Q9NQS3

Full Name:

poliovirus receptor-related 3

Calculated MW:

61 kDa

Purification Method:

Protein G purification

CloneNo.:

2C8C2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51308-1, Nectin 3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60892-1-PBS. Detection antibody: 60892-2-PBS. Standard:Ag28680. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP51308-2, Nectin 3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60892-2-PBS. Detection antibody: 60892-3-PBS. Standard:Ag28680. Range: 0.391-100 ng/mL