

For Research Use Only

# P130 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60893-2-PBS



## Basic Information

Catalog Number:

60893-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28262

GenBank Accession Number:

BC006769

GeneID (NCBI):

9221

UNIPROT ID:

Q14978

Full Name:

nucleolar and coiled-body

phosphoprotein 1

Calculated MW:

73 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1D8E6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP51309-1, P130 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60893-1-PBS. Detection antibody: 60893-2-PBS. Standard:Ag28262. Range: 0.781-200 ng/mL.



Cytometric bead array standard curve of MP51309-2, P130 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60893-3-PBS. Detection antibody: 60893-2-PBS. Standard:Ag28262. Range: 12.5-200 ng/mL.