

For Research Use Only

P130 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60893-3-PBS



Basic Information

Catalog Number:

60893-3-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28262

GenBank Accession Number:

BC006769

GeneID (NCBI):

9221

UNIPROT ID:

Q14978

Full Name:

nucleolar and coiled-body phosphoprotein 1

Calculated MW:

73 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2A10A5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51309-2, P130 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60893-3-PBS. Detection antibody: 60893-2-PBS. Standard:Ag28262. Range: 12.5-200 ng/mL