

For Research Use Only

TRIM31 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60895-2-PBS



Basic Information

Catalog Number:

60895-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29697

GenBank Accession Number:

BC017017

GeneID (NCBI):

11074

UNIPROT ID:

Q9BZY9

Full Name:

tripartite motif-containing 31

Calculated MW:

425 aa, 48 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2A11A2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51315-1, TRIM31 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60895-1-PBS. Detection antibody: 60895-2-PBS. Standard:Ag29697. Range: 1.563-200 ng/mL.