

For Research Use Only

# TRIM31 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60895-3-PBS



## Basic Information

Catalog Number:

60895-3-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG29697

GenBank Accession Number:

BC017017

GeneID (NCBI):

11074

UNIPROT ID:

Q9BZY9

Full Name:

tripartite motif-containing 31

Calculated MW:

425 aa, 48 kDa

Purification Method:

Protein A purification

CloneNo.:

1F3A3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51315-2, TRIM31 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60895-3-PBS. Detection antibody: 60895-4-PBS. Standard:Ag29697. Range: 1.563-200 ng/mL.