

For Research Use Only

BMP4 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60897-1-PBS



Basic Information

Catalog Number: 60897-1-PBS	GenBank Accession Number: BC020546	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 652	CloneNo.: 1B8E11
Source: Mouse	UNIPROT ID: P12644	
Isotype: IgG2b	Full Name: bone morphogenetic protein 4	
Immunogen Catalog Number: AG28840	Calculated MW: 47 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array of cell line supernatant using MP51307-1, BMP-4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60897-1-PBS. Detection antibody: 69015-1-PBS. Supernatant (50 μ L/well) from the following cell lines were assayed: HT-29, TF-1.



Cytometric bead array standard curve of MP51307-1, BMP-4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60897-1-PBS. Detection antibody: 69015-1-PBS. Standard:HZ-1045. Range: 6.25-100 ng/mL.