

For Research Use Only

EGF Monoclonal antibody, PBS Only

Catalog Number: 60900-1-PBS



Basic Information

Catalog Number:

60900-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG17529

GenBank Accession Number:

BC093731

GeneID (NCBI):

1950

UNIPROT ID:

P01133

Full Name:

epidermal growth factor (beta-urogastrone)

Calculated MW:

1207 aa, 134 kDa

Purification Method:

Protein A purification

CloneNo.:

3H9G11

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

