

For Research Use Only

ITLN1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60906-1-PBS



Basic Information

Catalog Number: 60906-1-PBS	GenBank Accession Number: BC020664	Purification Method: Protein A Magarose purification
Concentration: 1 mg/ml	GeneID (NCBI): 55600	CloneNo.: 1E9H1
Source: Mouse	UNIPROT ID: Q8WWAO	
Isotype: IgG3	Full Name: intelectin 1 (galactofuranose binding)	
Immunogen Catalog Number: AG28617	Calculated MW: 313 aa, 35 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51319-1, ITLN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60906-1-PBS. Detection antibody: 60906-2-PBS.



Cytometric bead array standard curve of MP51319-1, ITLN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60906-1-PBS. Detection antibody: 60906-2-PBS. Standard:Ag28617. Range: 0.098-100 ng/mL