

For Research Use Only

ITLN1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60906-2-PBS



Basic Information

Catalog Number:

60906-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28617

GenBank Accession Number:

BC020664

GeneID (NCBI):

55600

UNIPROT ID:

Q8WWAO

Full Name:

intelectin 1 (galactofuranose binding)

Calculated MW:

313 aa, 35 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3A3D7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Sample test of MP51319-1, ITLN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60906-1-PBS. Detection antibody: 60906-2-PBS.



Sample test of MP51319-2, ITLN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60906-3-PBS. Detection antibody: 60906-2-PBS.



Cytometric bead array standard curve of MP51319-1, ITLN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60906-1-PBS. Detection antibody: 60906-2-PBS. Standard:Ag28617. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP51319-2, ITLN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60906-3-PBS. Detection antibody: 60906-2-PBS. Standard:Ag28617. Range: 0.098-100 ng/mL