

For Research Use Only

JUN Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60907-1-PBS



Basic Information

Catalog Number:

60907-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC068522

GeneID (NCBI):

3725

UNIPROT ID:

P05412

Full Name:

jun oncogene

Calculated MW:

331 aa, 36 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1F12E3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51323-1, JUN Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66313-2-PBS. Detection antibody: 60907-1-PBS.



Cytometric bead array standard curve of MP51323-1, JUN Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66313-2-PBS. Detection antibody: 60907-1-PBS. Standard: Ag17419. Range: 3.906-500 ng/mL.