

For Research Use Only

EPHB2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60908-1-PBS



Basic Information

Catalog Number:	60908-1-PBS	GenBank Accession Number:	BC018763	Purification Method:	Protein G Magarose purification
Concentration:	1 mg/ml	GeneID (NCBI):	2048	CloneNo.:	1G2G1
Source:	Mouse	UNIPROT ID:	P29323		
Isotype:	IgG1	Full Name:	EPH receptor B2		
Immunogen Catalog Number:	AG28796	Calculated MW:	987 aa, 108 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51327-1, EPHB2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60908-1-PBS. Detection antibody: 60908-2-PBS.



Cytometric bead array standard curve of MP51327-1, EPHB2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60908-1-PBS. Detection antibody: 60908-2-PBS. Standard:Ag28796. Range: 7.813-500 ng/mL.