

For Research Use Only

MAPK11 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60913-3-PBS



Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
60913-3-PBS	BC027933	Protein A purification
Concentration:	GeneID (NCBI):	CloneNo.:
1 mg/ml	5600	2D2A6
Source:	UNIPROT ID:	
Mouse	Q15759	
Isotype:	Full Name:	
IgG2b	mitogen-activated protein kinase 11	
Immunogen Catalog Number:	Calculated MW:	
AG11986	364 aa, 41 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51335-2, MAPK11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60913-1-PBS. Detection antibody: 60913-3-PBS.



Cytometric bead array standard curve of MP51335-2, MAPK11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60913-1-PBS. Detection antibody: 60913-3-PBS. Standard:Ag11986. Range: 1.953-500 ng/mL