

For Research Use Only

CHAC2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60914-1-PBS



Basic Information

Catalog Number:

60914-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG9586

GenBank Accession Number:

BC019239

GeneID (NCBI):

494143

UNIPROT ID:

Q8WUX2

Full Name:

ChaC, cation transport regulator homolog 2 (E. coli)

Calculated MW:

184 aa, 21 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2C11E2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51336-1, CHAC2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60914-1-PBS. Detection antibody: 60914-2-PBS.



Cytometric bead array standard curve of MP51336-1, CHAC2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60914-1-PBS. Detection antibody: 60914-2-PBS. Standard:Ag9586. Range: 0.195-200 ng/mL