

For Research Use Only

# IGSF21 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60915-1-PBS



## Basic Information

Catalog Number: 60915-1-PBS	GenBank Accession Number: BC007618	Purification Method: Protein A Magarose purification
Concentration: 1 mg/ml	GeneID (NCBI): 84966	CloneNo.: 2A8C5
Source: Mouse	UNIPROT ID: Q96ID5	
Isotype: IgG2a	Full Name: immunoglobulin superfamily, member 21	
Immunogen Catalog Number: AG15062	Calculated MW: 467 aa, 52 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**


Selected Validation Data




Sample test of MP51337-1, IGSF21 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60915-1-PBS. Detection antibody: 60915-2-PBS.



Sample test of MP51337-2, IGSF21 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60915-1-PBS. Detection antibody: 60915-3-PBS.



Cytometric bead array standard curve of MP51337-1, IGSF21 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60915-1-PBS. Detection antibody: 60915-2-PBS. Standard:Ag15062. Range: 1.953-500 ng/mL



Cytometric bead array standard curve of MP51337-2, IGSF21 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60915-1-PBS. Detection antibody: 60915-3-PBS. Standard:Ag15062. Range: 1.953-500 ng/mL