

For Research Use Only

RBP3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60921-1-PBS



Basic Information

| | | |
|-------------------------------------|---|--|
| Catalog Number: 60921-1-PBS | GenBank Accession Number: BC039844 | Purification Method: Protein G purification |
| Concentration: 1 mg/ml | GeneID (NCBI): 5949 | CloneNo.: 1C8D10 |
| Source: Mouse | UNIPROT ID: P10745 | |
| Isotype: IgG1 | Full Name: retinol binding protein 3, interstitial | |
| Immunogen Catalog Number: AG6182 | Calculated MW: 135 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51339-1, RBP3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60921-1-PBS. Detection antibody: 60921-2-PBS.



Cytometric bead array standard curve of MP51339-1, RBP3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60921-1-PBS. Detection antibody: 60921-2-PBS. Standard: Ag6182. Range: 0.391-400 ng/mL.