

For Research Use Only

RBP3 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60921-2-PBS



Basic Information

Catalog Number: 60921-2-PBS	GenBank Accession Number: BC039844	Purification Method: Protein G Magarose purification
Concentration: 1 mg/ml	GeneID (NCBI): 5949	CloneNo.: 1A10C10
Source: Mouse	UNIPROT ID: P10745	
Isotype: IgG1	Full Name: retinol binding protein 3, interstitial	
Immunogen Catalog Number: AG6182	Calculated MW: 135 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51339-1, RBP3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60921-1-PBS. Detection antibody: 60921-2-PBS.



Cytometric bead array standard curve of MP51339-1, RBP3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60921-1-PBS. Detection antibody: 60921-2-PBS. Standard: Ag6182. Range: 0.391-400 ng/mL.