

For Research Use Only

FGF1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60927-1-PBS



Basic Information

Catalog Number: 60927-1-PBS	GenBank Accession Number: BC032697	Purification Method: Protein G purification
Concentration: 1 mg/ml	GeneID (NCBI): 2246	CloneNo.: 2C9A1
Source: Mouse	UNIPROT ID: P05230	
Isotype: IgG1	Full Name: fibroblast growth factor 1 (acidic)	
Immunogen Catalog Number: AG11672	Calculated MW: 155 aa, 17 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51346-1, FGF1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60927-1-PBS. Detection antibody: 60927-2-PBS.



Cytometric bead array standard curve of MP51346-1, FGF1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60927-1-PBS. Detection antibody: 60927-2-PBS. Standard: Ag11672. Range: 0.391-400 ng/mL