

For Research Use Only

# GLP1R Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60931-2-PBS



## Basic Information

Catalog Number: 60931-2-PBS	GenBank Accession Number: BC112126	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 2740	CloneNo.: 1A7C3
Source: Mouse	UNIPROT ID: P43220	
Isotype: IgG2a	Full Name: glucagon-like peptide 1 receptor	
Immunogen Catalog Number: AG23529	Calculated MW: 463 aa, 53 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA, Sample test

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

Selected Validation Data



Sample test of MP51365-1, GLP1R Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60931-1-PBS. Detection antibody: 60931-2-PBS.



Cytometric bead array standard curve of MP51365-1, GLP1R Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60931-1-PBS. Detection antibody: 60931-2-PBS. Standard:Ag23529. Range: 0.391-200 ng/mL