

For Research Use Only

S100A7 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60932-2-PBS



Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
60932-2-PBS	BC034687	Protein G Magarose purification
Concentration:	GeneID (NCBI):	CloneNo.:
1 mg/ml	6278	1C1A5
Source:	UNIPROT ID:	
Mouse	P31151	
Isotype:	Full Name:	
IgG1	S100 calcium binding protein A7	
Immunogen Catalog Number:	Calculated MW:	
AG23611	101 aa, 11 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51348-1, S100A7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60932-1-PBS. Detection antibody: 60932-2-PBS.



Cytometric bead array standard curve of MP51348-1, S100A7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60932-1-PBS. Detection antibody: 60932-2-PBS. Standard:Ag23611. Range: 0.391-400 ng/mL.