

For Research Use Only

SMARCAL1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60933-1-PBS



Basic Information

Catalog Number:

60933-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG23707

GenBank Accession Number:

BC016482

GeneID (NCBI):

50485

UNIPROT ID:

Q9NZC9

Full Name:

SWI/SNF related, matrix associated,
actin dependent regulator of
chromatin, subfamily a-like 1

Calculated MW:

954 aa, 106 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3H2F4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51350-1, SMARCAL1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60933-1-PBS. Detection antibody: 60933-2-PBS.



Cytometric bead array standard curve of MP51350-1, SMARCAL1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60933-1-PBS. Detection antibody: 60933-2-PBS. Standard:Ag23707. Range: 0.098-100 ng/mL.