

For Research Use Only

NRAS Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60934-1-PBS



Basic Information

Catalog Number: 60934-1-PBS	GenBank Accession Number: BC005219	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 4893	CloneNo.: 2B11F3
Source: Mouse	UNIPROT ID: P01111	
Isotype: IgG2a	Full Name: neuroblastoma RAS viral (v-ras) oncogene homolog	
Immunogen Catalog Number: AG29746	Calculated MW: 21 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51352-1, NRAS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60934-1-PBS. Detection antibody: 60934-2-PBS.



Cytometric bead array standard curve of MP51352-1, NRAS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60934-1-PBS. Detection antibody: 60934-2-PBS. Standard:Ag29746. Range: 0.391-400 ng/mL.