

For Research Use Only

# NRAS Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60934-2-PBS



## Basic Information

Catalog Number:

60934-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG29746

GenBank Accession Number:

BC005219

GeneID (NCBI):

4893

UNIPROT ID:

P01111

Full Name:

neuroblastoma RAS viral (v-ras)  
oncogene homolog

Calculated MW:

21 kDa

Purification Method:

Protein A purification

CloneNo.:

1A6C10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Sample test of MP51352-1, NRAS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60934-1-PBS. Detection antibody: 60934-2-PBS.



Cytometric bead array standard curve of MP51352-1, NRAS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60934-1-PBS. Detection antibody: 60934-2-PBS. Standard:Ag29746. Range: 0.391-400 ng/mL.