

For Research Use Only

NOV Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60935-1-PBS



Basic Information

Catalog Number: 60935-1-PBS	GenBank Accession Number: BC015028	Purification Method: Protein G purification
Concentration: 1 mg/ml	GeneID (NCBI): 4856	CloneNo.: 2H8D4
Source: Mouse	UNIPROT ID: P48745	
Isotype: IgG1	Full Name: nephroblastoma overexpressed gene	
Immunogen Catalog Number: AG30138	Calculated MW: 211 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array of Human urine using MP51353-1, NOV Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60935-1-PBS. Detection antibody: 60935-2-PBS.



Sample test of MP51353-1, NOV Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60935-1-PBS. Detection antibody: 60935-2-PBS.



Cytometric bead array standard curve of MP51353-1, NOV Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60935-1-PBS. Detection antibody: 60935-2-PBS. Standard: Ag30138. Range: 0.098-100 ng/mL.