

For Research Use Only

AKAP12 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60936-1-PBS



Basic Information

Catalog Number: 60936-1-PBS	GenBank Accession Number: BC157067	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 9590	CloneNo.: 2C2D11
Source: Mouse	UNIPROT ID: Q02952	
Isotype: IgG2a	Full Name: A kinase (PRKA) anchor protein 12	
Immunogen Catalog Number: AG32851	Calculated MW: 1782 aa, 191 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51354-1, AKAP12 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60936-1-PBS. Detection antibody: 60936-2-PBS.



Cytometric bead array standard curve of MP51354-1, AKAP12 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60936-1-PBS. Detection antibody: 60936-2-PBS. Standard:Ag32851. Range: 0.098-100 ng/mL.