

For Research Use Only

# Aiolos/IKZF3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60938-1-PBS



## Basic Information

Catalog Number:

60938-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5537

GenBank Accession Number:

BC032707

GeneID (NCBI):

22806

UNIPROT ID:

Q9UKT9

Full Name:

IKAROS family zinc finger 3 (Aiolos)

Calculated MW:

509 aa, 58 kDa

Purification Method:

Protein G purification

CloneNo.:

3D7E2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Sample test of MP51356-1, Aiolos/IKZF3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60938-1-PBS. Detection antibody: 60938-2-PBS.



Cytometric bead array standard curve of MP51356-1, Aiolos/IKZF3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60938-1-PBS. Detection antibody: 60938-2-PBS. Standard:Ag5537. Range: 0.391-400 ng/mL.