

For Research Use Only

# MCM3 Monoclonal antibody, PBS Only

Catalog Number: 60939-1-PBS



## Basic Information

Catalog Number:

60939-1-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG8122

GenBank Accession Number:

BC001626

GeneID (NCBI):

4172

UNIPROT ID:

P25205

Full Name:

minichromosome maintenance  
complex component 3

Calculated MW:

91 kDa

Observed MW:

102-105 kDa

Purification Method:

Protein G purification

CloneNo.:

2G9C5

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

## Background Information

MCM3, also named as HCC5, P1-MCM3, p102, RLF subunit beta and RLFB, belongs to the MCM family. MCM3 is essential for the initiation of DNA replication and is involved in ensuring that DNA replication is initiated precisely once per cell cycle(PMID: 33093917). The overexpression of MCM3 is related to proliferation of tumor cells(PMID: 31186514).

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 60939-1-Ig (MCM3 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 60939-1-PBS in a different storage buffer formulation.