

For Research Use Only

# GMFG Monoclonal antibody, PBS Only

Catalog Number: 60941-1-PBS



## Basic Information

Catalog Number: 60941-1-PBS	GenBank Accession Number: BC032819	Purification Method: Protein G purification
Concentration: 1 mg/ml	GeneID (NCBI): 9535	CloneNo.: 2B3E2
Source: Mouse	UNIPROT ID: O60234	
Isotype: IgG1	Full Name: glia maturation factor, gamma	
Immunogen Catalog Number: AG4538	Calculated MW: 141 aa, 17 kDa	

## Applications

Tested Applications:  
Indirect ELISA  
Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
Storage Buffer:  
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Selected Validation Data