

For Research Use Only

GMFG Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60941-3-PBS



Basic Information

Catalog Number: 60941-3-PBS	GenBank Accession Number: BC032819	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 9535	CloneNo.: 1D12C8
Source: Mouse	UNIPROT ID: O60234	
Isotype: IgG3	Full Name: glia maturation factor, gamma	
Immunogen Catalog Number: AG4538	Calculated MW: 141 aa, 17 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51410-1, GMFG Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60941-2-PBS. Detection antibody: 60941-3-PBS.



Cytometric bead array standard curve of MP51410-1, GMFG Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60941-2-PBS. Detection antibody: 60941-3-PBS. Standard: Ag4538. Range: 0.781-100 ng/mL.