

For Research Use Only

PLAT Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60942-2-PBS



Basic Information

Catalog Number: 60942-2-PBS	GenBank Accession Number: BC095403	Purification Method: Protein G Magarose purification
Source: Mouse	GeneID (NCBI): 5327	CloneNo.: 1E12E9
Isotype: IgG1	UNIPROT ID: P00750	
Immunogen Catalog Number: EG30025	Full Name: plasminogen activator, tissue	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51573-1, PLAT Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60942-1-PBS. Detection antibody: 60942-2-PBS.



Cytometric bead array standard curve of MP51573-1, PLAT Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60942-1-PBS. Detection antibody: 60942-2-PBS. Standard: Eg30025. Range: 0.781-100 ng/mL.