

For Research Use Only

IGSF8 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60946-2-PBS



Basic Information

Catalog Number:

60946-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5994

GenBank Accession Number:

BC053881

GeneID (NCBI):

93185

UNIPROT ID:

Q969P0

Full Name:

immunoglobulin superfamily,

member 8

Calculated MW:

65 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2F11G2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com


This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array of Human urine using MP51358-1, IGSF8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60946-1-PBS. Detection antibody: 60946-2-PBS.

Sample test of MP51358-1, IGSF8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60946-1-PBS. Detection antibody: 60946-2-PBS.



Cytometric bead array standard curve of MP51358-1, IGSF8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60946-1-PBS. Detection antibody: 60946-2-PBS. Standard:Ag5994. Range: 0.781-100 ng/mL