

For Research Use Only

SERTAD1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60947-1-PBS



Basic Information

Catalog Number: 60947-1-PBS	GenBank Accession Number: BC002670	Purification Method: Protein G purification
Concentration: 1 mg/ml	GeneID (NCBI): 29950	CloneNo.: 1E6D6
Source: Mouse	UNIPROT ID: Q9UHV2	
Isotype: IgG1	Full Name: SERTA domain containing 1	
Immunogen Catalog Number: AG25568	Calculated MW: 25 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51359-1, SERTAD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60947-1-PBS. Detection antibody: 60947-2-PBS.

Cytometric bead array standard curve of MP51359-1, SERTAD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60947-1-PBS. Detection antibody: 60947-2-PBS. Standard: Ag25568. Range: 0.098-100 ng/mL