

For Research Use Only

# SERTAD1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60947-2-PBS



## Basic Information

Catalog Number: 60947-2-PBS	GenBank Accession Number: BC002670	Purification Method: Protein G Magarose purification
Concentration: 1 mg/ml	GeneID (NCBI): 29950	CloneNo.: 3E11A9
Source: Mouse	UNIPROT ID: Q9UHV2	
Isotype: IgG1	Full Name: SERTA domain containing 1	
Immunogen Catalog Number: AG25568	Calculated MW: 25 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA, Sample test

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:



T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sample test of MP51359-1, SERTAD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60947-1-PBS. Detection antibody: 60947-2-PBS.

Cytometric bead array standard curve of MP51359-1, SERTAD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60947-1-PBS. Detection antibody: 60947-2-PBS. Standard: Ag25568. Range: 0.098-100 ng/mL