

For Research Use Only

# PTPRE Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60950-2-PBS



## Basic Information

<b>Catalog Number:</b> 60950-2-PBS	<b>GenBank Accession Number:</b> BC050062	<b>Purification Method:</b> Protein G Magarose purification
<b>Concentration:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 5791	<b>CloneNo.:</b> 1F5C11
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> P23469	
<b>Isotype:</b> IgG1	<b>Full Name:</b> protein tyrosine phosphatase, receptor type, E	
<b>Immunogen Catalog Number:</b> AG4932	<b>Calculated MW:</b> 81 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA, Sample test

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sample test of MP51362-1, PTPRE Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60950-1-PBS. Detection antibody: 60950-2-PBS.



Cytometric bead array standard curve of MP51362-1, PTPRE Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60950-1-PBS. Detection antibody: 60950-2-PBS. Standard: Ag4932. Range: 0.195-200 ng/mL