

For Research Use Only

HSD17B4 Monoclonal antibody, PBS Only (Capture)

Catalog Number:60951-1-PBS



Basic Information

Catalog Number: 60951-1-PBS	GenBank Accession Number: BC003098	Purification Method: Protein A purification
Source: Mouse	GeneID (NCBI): 3295	CloneNo.: 1A3H11
Isotype: IgG2a	UNIPROT ID: P51659	
Immunogen Catalog Number: AG7221	Full Name: hydroxysteroid (17-beta) dehydrogenase 4	
	Calculated MW: 80 kDa	
	Observed MW: 80 kDa, 45 kDa	

Applications

Tested Applications:
WB, IF/ICC, Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

HSD17B4 (17-beta-hydroxysteroid dehydrogenase 4) is also named as Peroxisomal multifunctional enzyme type 2, D-bifunctional protein or multifunctional protein 2. It codes for a 80 kDa enzyme containing three distinct functional domains and is localized in peroxisomes. It is a bifunctional enzyme acting on the peroxisomal beta-oxidation pathway for fatty acids and catalyzing the formation of 3-ketoacyl-CoA intermediates from both straight-chain and 2-methyl-branched-chain fatty acids. After peroxisomal import, the full-length protein is proteolytically cleaved to yield a 35-kDa dehydrogenase subunit and a 45-kDa hydratase subunit containing the hydratase and SCP domains (PMID: 28868548, 24602372).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 60951-1-Ig (HSD17B4 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 60951-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using HSD17B4 antibody (60951-1-Ig, Clone: 1A3H11) at dilution of 1:2000 and CoraLite@488-Conjugated Goat Anti-Mouse IgG(H+L) (SA00013-1), CL594-phalloidin (red). This data was developed using the same antibody clone with 60951-1-PBS in a different storage buffer formulation.



Sample test of MP51363-1, HSD17B4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60951-1-PBS. Detection antibody: 60951-2-PBS.



Cytometric bead array standard curve of MP51363-1, HSD17B4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60951-1-PBS. Detection antibody: 60951-2-PBS. Standard: Ag7221. Range: 0.391-200 ng/mL.