

For Research Use Only

GFM2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60952-1-PBS



Basic Information

Catalog Number:

60952-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG36976

GenBank Accession Number:

BC015712

GeneID (NCBI):

84340

UNIPROT ID:

Q96959

Full Name:

G elongation factor, mitochondrial 2

Calculated MW:

779 aa, 87 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1G1F8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51366-1, GFM2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60952-1-PBS. Detection antibody: 60952-2-PBS.



Cytometric bead array standard curve of MP51366-1, GFM2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60952-1-PBS. Detection antibody: 60952-2-PBS. Standard: Ag36976. Range: 12.5-400 ng/mL.