

For Research Use Only

LCORL Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60954-1-PBS



Basic Information

Catalog Number:

60954-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG11344

GenBank Accession Number:

BC037322

GeneID (NCBI):

254251

UNIPROT ID:

Q8N3X6

Full Name:

ligand dependent nuclear receptor
corepressor-like

Calculated MW:

318 aa, 35 kDa, 67 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3D5E11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51368-1, LCORL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60954-1-PBS. Detection antibody: 60954-2-PBS.

Cytometric bead array standard curve of MP51368-1, LCORL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60954-1-PBS. Detection antibody: 60954-2-PBS. Standard: Ag11344. Range: 0.195-50 ng/mL.