

For Research Use Only

# LCORL Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60954-2-PBS



## Basic Information

Catalog Number:

60954-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG11344

GenBank Accession Number:

BC037322

GeneID (NCBI):

254251

UNIPROT ID:

Q8N3X6

Full Name:

ligand dependent nuclear receptor  
corepressor-like

Calculated MW:

318 aa, 35 kDa, 67 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3A11F2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:



T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sample test of MP51368-1, LCORL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60954-1-PBS. Detection antibody: 60954-2-PBS.

Cytometric bead array standard curve of MP51368-1, LCORL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60954-1-PBS. Detection antibody: 60954-2-PBS. Standard: Ag11344. Range: 0.195-50 ng/mL.