

For Research Use Only

HDAC7 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60956-2-PBS



Basic Information

| | | |
|--|--|---|
| Catalog Number: 60956-2-PBS | GenBank Accession Number: BC006453 | Purification Method: Protein A purification |
| Concentration: 1 mg/ml | GeneID (NCBI): 51564 | CloneNo.: 1C4G5 |
| Source: Mouse | UNIPROT ID: Q8WUI4 | |
| Isotype: IgG2b | Full Name: histone deacetylase 7 | |
| Immunogen Catalog Number: AG8277 | Calculated MW: 103 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51370-1, HDAC7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60956-1-PBS. Detection antibody: 60956-2-PBS.



Cytometric bead array standard curve of MP51370-1, HDAC7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60956-1-PBS. Detection antibody: 60956-2-PBS. Standard: Ag8277. Range: 0.098-100 ng/mL