

For Research Use Only

P19 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60958-1-PBS



Basic Information

Catalog Number:

60958-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28546

GenBank Accession Number:

BC001822

GeneID (NCBI):

1032

UNIPROT ID:

P55273

Full Name:

cyclin-dependent kinase inhibitor 2D
(p19, inhibits CDK4)

Calculated MW:

19 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1F5H10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51373-1, P19 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60958-1-PBS. Detection antibody: 60958-2-PBS.

Cytometric bead array standard curve of MP51373-1, P19 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60958-1-PBS. Detection antibody: 60958-2-PBS. Standard: Ag28546. Range: 0.391-200 ng/mL.