

For Research Use Only

TAF11 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60962-2-PBS



Basic Information

Catalog Number: 60962-2-PBS	GenBank Accession Number: BC021972	Purification Method: Protein A Magarose purification
Concentration: 1 mg/ml	GeneID (NCBI): 6882	CloneNo.: 3B10G9
Source: Mouse	UNIPROT ID: Q15544	
Isotype: IgG2a	Full Name: TAF11 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 28kDa	
Immunogen Catalog Number: AG9629	Calculated MW: 211 aa, 23 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51378-1, TAF11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60962-1-PBS. Detection antibody: 60962-2-PBS.

Cytometric bead array standard curve of MP51378-1, TAF11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60962-1-PBS. Detection antibody: 60962-2-PBS. Standard: Ag9629. Range: 12.5-200 ng/mL.