

For Research Use Only

NISCH Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60964-2-PBS



Basic Information

Catalog Number: 60964-2-PBS	GenBank Accession Number: BC038102	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 11188	CloneNo.: 3D8E2
Source: Mouse	UNIPROT ID: Q9Y2I1	
Isotype: IgG2a	Full Name: nischarin	
Immunogen Catalog Number: AG33989	Calculated MW: 1504 aa, 167 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51380-1, NISCH Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60964-1-PBS. Detection antibody: 60964-2-PBS.



Cytometric bead array standard curve of MP51380-1, NISCH Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60964-1-PBS. Detection antibody: 60964-2-PBS. Standard: Ag33989. Range: 1.563-400 ng/mL.