

For Research Use Only

# DDHD2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60965-1-PBS



## Basic Information

Catalog Number: 60965-1-PBS	GenBank Accession Number: BC010504	Purification Method: Protein G Magarose purification
Concentration: 1 mg/ml	GeneID (NCBI): 23259	CloneNo.: 3G3F8
Source: Mouse	UNIPROT ID: O94830	
Isotype: IgG1	Full Name: DDHD domain containing 2	
Immunogen Catalog Number: AG18391	Calculated MW: 711 aa, 81 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA, Sample test

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sample test of MP51382-1, DDHD2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60965-1-PBS. Detection antibody: 60965-2-PBS.



Cytometric bead array standard curve of MP51382-1, DDHD2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60965-1-PBS. Detection antibody: 60965-2-PBS. Standard: Ag18391. Range: 3.125-400 ng/mL