

For Research Use Only

FHL3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60966-1-PBS



Basic Information

Catalog Number: 60966-1-PBS	GenBank Accession Number: BC011697	Purification Method: Protein A purification
Source: Mouse	GeneID (NCBI): 2275	CloneNo.: 1G3A10
Isotype: IgG2b	UNIPROT ID: Q13643	
Immunogen Catalog Number: AG1495	Full Name: four and a half LIM domains 3	
	Calculated MW: 31 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Sample test of MP51383-1, FHL3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60966-1-PBS. Detection antibody: 60966-2-PBS.



Cytometric bead array standard curve of MP51383-1, FHL3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60966-1-PBS. Detection antibody: 60966-2-PBS. Standard: Ag1495. Range: 6.25-400 ng/mL.



Biolayer interferometry (BLI) kinetic assays of 60966-1-Ig against Human FHL3 were performed. The affinity constant is 0.78 nM.