

For Research Use Only

FHL3 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60966-2-PBS



Basic Information

Catalog Number:

60966-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG1495

GenBank Accession Number:

BC011697

GeneID (NCBI):

2275

UNIPROT ID:

Q13643

Full Name:

four and a half LIM domains 3

Calculated MW:

31 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3D7A12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51383-1, FHL3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60966-1-PBS. Detection antibody: 60966-2-PBS.



Cytometric bead array standard curve of MP51383-1, FHL3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60966-1-PBS. Detection antibody: 60966-2-PBS. Standard:Ag1495. Range: 6.25-400 ng/mL.