

For Research Use Only

TXLNA Monoclonal antibody, PBS Only

Catalog Number:60968-1-PBS



Basic Information

Catalog Number: 60968-1-PBS	GenBank Accession Number: BC103823	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 200081	CloneNo.: 1D7G6
Isotype: IgG1	UNIPROT ID: P40222	
Immunogen Catalog Number: AG25020	Full Name: taxilin alpha	
	Calculated MW: 546 aa, 62 kDa	

Applications

Tested Applications:
Indirect ELISA
Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data