

For Research Use Only

TMPRSS4 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60969-1-PBS



Basic Information

Catalog Number:

60969-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG1804

GenBank Accession Number:

BC011703

GeneID (NCBI):

56649

UNIPROT ID:

Q9NRS4

Full Name:

transmembrane protease, serine 4

Calculated MW:

48 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2E10E8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51391-1, TMPRSS4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60969-1-PBS. Detection antibody: 60969-2-PBS.



Cytometric bead array standard curve of MP51391-1, TMPRSS4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60969-1-PBS. Detection antibody: 60969-2-PBS. Standard:Ag1804. Range: 0.781-100 ng/mL.