

For Research Use Only

SYNJ1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60973-2-PBS



Basic Information

Catalog Number: 60973-2-PBS	GenBank Accession Number: BC136603	Purification Method: Protein G Magarose purification
Concentration: 1 mg/ml	GeneID (NCBI): 8867	CloneNo.: 1C8B1
Source: Mouse	UNIPROT ID: O43426	
Isotype: IgG1	Full Name: synaptojanin 1	
Immunogen Catalog Number: AG17719	Calculated MW: 1612 aa, 178 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51398-1, SYNJ1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60973-1-PBS. Detection antibody: 60973-2-PBS.

Cytometric bead array standard curve of MP51398-1, SYNJ1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60973-1-PBS. Detection antibody: 60973-2-PBS. Standard: Ag17719. Range: 0.195-50 ng/mL.