

For Research Use Only

TBC1D17 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60975-1-PBS



Basic Information

Catalog Number: 60975-1-PBS	GenBank Accession Number: BC003516	Purification Method: Protein G purification
Concentration: 1 mg/ml	GeneID (NCBI): 79735	CloneNo.: 1G9G1
Source: Mouse	UNIPROT ID: Q9HA65	
Isotype: IgG1	Full Name: TBC1 domain family, member 17	
Immunogen Catalog Number: AG14358	Calculated MW: 648 aa, 73 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51406-1, TBC1D17 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60975-1-PBS. Detection antibody: 60975-2-PBS.

Cytometric bead array standard curve of MP51406-1, TBC1D17 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60975-1-PBS. Detection antibody: 60975-2-PBS. Standard: Ag14358. Range: 1.563-200 ng/mL